Professional Veterinary Diagnostic Test
BioPronix product line
 For Veterinary use only
 For in vitro use only
 Manufactured by Agrolabo S.p.A.

INTERPRETATION OF THE RESULTS

Negative Result

The test is negative if a red line appears in window C (control line).

Positive Result Ehrlichia

The test is positive for Ehrlichia if two red lines appear in column 1. One line band will appear in window T (test line) and one in window C (control line).

Positive Result Anaplasma

The test is positive for Anaplasma if two red lines appear in column 2. One line band will appear in window T (test line) and one in window C (control line).

Invalid Result

The test is invalid if no coloured band appears in window C (control line), even if a coloured band appears in window T (test line).

WARNING

Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

QUICK GUIDE



1. Remove the device from its wrapping and place it on a flat surface.

Use the pipette to add 1 drop of serum, plasma or whole blood with anticoagulant in both windows S of the device

2. Hold the diluent bottle in a vertical position and add 2 drops of diluent in both windows S of the device

3. Wait 10 minutes to read the results, not later than 15 minutes

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

EHRlichia/ANAPLASMA IC - Rapid immunochromatographic test for the detection of antibodies to Ehrlichia canis and Anaplasma spp.

TEST PRINCIPLE

Ehrlichia/Anaplasma IC is an immunochromatographic test for the detection of antibodies to Ehrlichia canis and Anaplasma spp. (A. phagocytophilum / A. platys) in canine whole blood, serum or plasma.

If the sample contains specific antibodies against Ehrlichia canis and Anaplasma spp. (positive sample), these antibodies will bind to the colloidal gold conjugated Ehrlichia and Anaplasma antigens included in the pad of the test, forming the Ag-Ab complex. The complex Ag-Ab migrates along the membrane and will bind again to the Ehrlichia or Anaplasma antigen at T window level forming a red line (test line).

The sample continues to migrate along the membrane up to window C, where a reagent able to bind the antigen conjugated with colloidal gold is immobilized resulting in the formation of a red line (control line).

KIT COMPONENTS

- Disposable pipettes
- Vial with diluent
- EHRlichia/ANAPLASMA IC devices

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C).

When properly stored, the product remains stable for its full shelf life.

Keep the kit away from excessive heat/cold or direct sunlight.

PRECAUTIONS

IMPORTANT: open individual device packaging immediately prior to use.

Do not use kit components after expiry date.

The test device should be placed horizontally on a flat surface while testing.

Use a new disposable pipette for each sample.

SAMPLE

All samples need to be brought to room temperature before use. Previously frozen or stored samples must be centrifuged before use.

Whole blood samples with anticoagulant (EDTA, heparin, etc) may be used. Whole blood samples should be fresh or stored in a refrigerator for no longer than 5 days.

Serum or plasma that is fresh, or that has been refrigerated (+2-8°C) for no longer than 5 days, can be used. For longer storage, samples should be frozen.

Haemolyzed samples may affect test results.

TEST PROCEDURE

1. Remove the device from its individual wrapping and place it on a flat surface.

2. Gently mix the sample by inverting. Holding the pipette approximately 1 cm from the device, add 1 drop of serum, plasma or whole blood with anticoagulant in both windows S. Wait for the sample to be absorbed.

3. Holding the diluent bottle in a vertical position, promptly add 2 drops of the diluent buffer in both windows S.

4. Read the result after 10 minutes and no later than 15 minutes.