FLUO ENCEPHALITOZOON cuniculi

Kit IFA for the detection of anti-Encephalitozoon cuniculi IgG antibodies in rabbit serum or plasma

INSTRUCTIONS FOR USE

INTENDED USE

The FLUO ENCEPHALITOZOON cuniculi kit is intended for the semi-quantitative detection of IgG antibodies to Encephalitozoon cuniculi in rabbit serum or plasma by indirect immunofluorescence assay.

TEST PRINCIPLE

Rabbit sera are diluted in buffered saline and incubated in the individual slide wells to allow reaction of patient antibody with Encephalitozoon cuniculi antigens fixed on the wells. Slides are then washed to remove unreacted serum proteins, a fluorescence-labelled FITC anti-rabbit IgG (conjugate) is added. This conjugate is allowed time to react with antigen-antibody complexes. The slides are washed again to remove unreacted conjugate. The resulting reactions can be visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescence of the Encephalitozoon cuniculi spores. In a negative reaction the apple-green fluorescence is not present or there is a fluorescent pattern different from the positive control well. Positive reactions may then be retested at a higher dilution to determine the highest reactive (endpoint titer).

TEST COMPONENTS

- 10 x 10 wells substrate slides;
- 1 vial of FITC anti-rabbit IgG conjugate-ready for use;
- 1 vial of positive control-ready for use;
- 1 vial of negative control- ready for use;
- 1 vial of mounting medium-ready for use;

MATERIALS REQUIRED BUT NOT SUPPLIED

PBS buffer (phosphate buffered saline 1X pH 7.2-7.4); test tubes for serum dilutions; precision pipettes and tips; 24x50 mm coverslips; fluorescence microscope with filter system for FITC (fluorescein isothiocyanate, excitation wavelength 465-495, barrier filter 515-555) and 400X magnification; +37°C incubator; humid chamber for slides.

STORAGE

Kit components should be stored at +2-8°C. Bring them to room temperature (+20-25 °C) before performing the test.

SAMPLE

Rabbit serum or plasma.


Allow sample to clot and separate serum by centrifugation. Transfer serum aseptically to a tightly closing sterile container. Store at +2-8 °C. If testing is to be delayed longer than 5 days, freezing the sample at -20 °C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two- and four-weeks intervals to check for titer changes.

TEST PROCEDURE

- Prepare 1:80 screening dilution for all untested sera. For sera found positive on a previous assay run, prepare 1:2 serial dilutions in PBS starting with 1:80.
- The positive and negative controls are ready to use. Do not dilute them. For each serum dilution to be tested add 20 µl to one slide well and record the location for later reference. For each assay run include the negative control and the positive control (20 µl).
- Place slides in a humid chamber and incubate for 30 minutes at +37°C.
- Remove humid chamber from incubator.
- Washing step: Tap remaining serum dilutions gently from the slides and shake the slides gently for 5 minutes in PBS. Repeat this step for another 5 minutes with fresh PBS. Briefly rinse the slides with distilled water. Tap remaining water gently from the slides and, if necessary, dry the Teflon mask between the wells with absorbent paper or cotton bud sticks. However, do not allow the antigen wells to dry out. If using a washing bottle, do not focus the stream directly onto the antigen wells.
- Add to each slide well 1 drop (20 µl) of conjugate, then return slides to the humid chamber for another 30 minutes incubation at +37°C. Incubation should be in the dark to protect the photosensitive conjugate.
- Repeat washing steps as previously described.
- Add 2 drops of mounting fluid to each slide and place the coverslip on, carefully removing air bubbles caught under the coverslip.
- Evaluate the slides at 400X magnification, comparing each well to the visual intensity and appearance of the patterns seen in the positive and negative control wells. Slides may be stored at +2-8°C in the dark for up to 7 days.

RESULTS

For the evaluation, a fluorescence microscope with a filter system for FITC and 400X magnification is required. The fluorescence pattern (form, density etc.) of the negative and positive control is considered as reference pattern. Other patterns of reactivity different than those seen in the controls must be considered non-specific, which means negative.

Negative test result at 1:80: report as negative for Encephalitozoon cuniculi IgG antibodies. The Encephalitozoon spores do not show a clear apple-green fluorescence.

Positive test result at dilution \geq 1:80: the Encephalitozoon spores show a bright, sharp and clear yellow-green fluorescence.

A further dilution of the positive samples is recommended to determine the highest dilution that is still positive (endpoint titer). An acute infection (2 to 4-fold titer increase) can only be determined by the titer determination of a coupled serum test (2 samples in an interval of 2-3 weeks).

The prevalence of Encephalitozoon cuniculi antibodies varies depending upon the geographic region and population being tested. Serological test results should be always interpreted in conjunction with case history and clinical symptoms. **Possible cross reactions:** due to cross reactivity between different types of Encephalitozoon cuniculi, it is not possible with IFA technique to distinguish between antibodies specific for Encephalitozoon cuniculi type I, II and/or III. Dog and cat: Encephalitozoon cuniculi single slides can be used also for the detection of antibodies anti-Encephalitozoon type II and III in canine and feline samples (using specific anti-specie FITC conjugates).


QUALITY CONTROL

The negative control serum and the positive control serum should be assayed with each daily run. The negative control well is an example of a non-reactive serum, with either uniform red counterstain or slight, but uniform greenish staining. The fluorescence intensity of the positive control is an example of medium intensity fluorescence. If one of the controls does not react as specified, the assay run should be considered invalid, reagent components and procedural steps should be re-checked, and the assay repeated.

PRECAUTIONS

For veterinary use only. Follow the instructions carefully. Since no testing can assure the absence of infectious agents, these reagents, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly. Conjugate is photosensitive, for protection store in the dark and return to storage after use. Conjugate contains Evans blue dye, avoid ingestion and skin contact.

Do not use after the expiration date.

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