FLUO LEISHMANIA

Kit IFA for the detection of anti-Leishmania IgG antibodies
in canine serum or plasma

INSTRUCTIONS FOR USE

INTENDED USE

The FLUO LEISHMANIA kit is intended for the semiquantitative detection of IgG antibodies to Leishmania in canine serum or plasma by indirect immunofluorescence assay

TEST PRINCIPLE

Canine sera are diluted in buffered saline and incubated in the individual slide wells to allow reaction of patient antibody with Leishmania coated on the wells. Slides are then washed to remove unreacted serum proteins, a fluorescence-labelled FITC anti-canine IgG (conjugate) is added. This conjugate is allowed time to react with antigen-antibody complexes. The slides are washed again to remove unreacted conjugate. The resulting reactions can be visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescent Leishmania membrane in each field. A negative reaction is seen with fluorescence unlike that seen in the positive control well. Positive reactions may then be retested at a higher dilution to determine the highest reactive (endpoint titer).

TEST COMPONENTS

- 10 x 10 wells substrate slides;
- 1 vial of FITC anti-dog IgG conjugate-ready for use;
- 1 vial of positive control-ready for use;
- 1 vial of negative control- ready for use;
- 1 vial of mounting medium-ready for use;

MATERIALS REQUIRED BUT NOT SUPPLIED

PBS buffer (phosphate buffered saline 1X pH 7,2-7,4); test tubes for serum dilutions; precision pipettes and tips; 24x50 mm coverslips; fluorescence microscope with filter system for FITC (fluorescein isothiocyanate, excitation wavelength 465-495, barrier filter 515-555) and 400X magnification; +37°C incubator; humid chamber for slides.

STORAGE

Kit components should be stored at +2-8°C. Bring them to room temperature (+20-25 °C) before performing the test.

SAMPLE

Canine serum or plasma.

Allow sample to clot and separate serum by centrifugation. Transfer serum aseptically to a tightly closing sterile container.

Store at +2-8°C. If testing is to be delayed longer than 5 days, freezing the sample at -20°C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two and four weeks intervals to check for titer changes.


TEST PROCEDURE

- Prepare **1:80** screening dilution for all untested sera. For sera found positive on a previous assay run, prepare serial two-fold dilution in PBS.
- The positive and negative controls are ready to use. Do not dilute them. For each serum dilution to be tested add 20 µl to one slide well and record the location for later reference. For each assay run include the negative control and the positive control (20 µl).
- Place slides in a humid chamber and incubate for 30 minutes at +37 °C.
- Remove humid chamber from incubator.
- Washing step: Tap remaining serum dilutions gently from the slides and shake the slides gently for 5 minutes in PBS. Repeat this step for another 5 minutes with fresh PBS. Briefly rinse the slides with distilled water. Tap remaining water gently from the slides and, if necessary, dry the Teflon mask between the wells with absorbent paper or cotton bud sticks. However, do not allow the antigen wells to dry out. If using a washing bottle, do not focus the stream directly onto the antigen wells.
- Add to each slide well 1 drop (20 µl) of conjugate, then return slides to the humid chamber for another 30 minutes incubation at +37°C. Incubation should be in the dark to protect the photosensitive conjugate.
- Repeat washing steps as previously described.
- Add 2 drops of mounting fluid to each slide and place the coverslip on, carefully removing air bubbles caught under the coverslip.
- Evaluate the slides at 400X magnification, comparing each well to the visual intensity and appearance of the Leishmania seen in the positive and negative control wells. Slides may be stored at +2-8°C in the dark for up to 7 days.

RESULTS

For the evaluation, a fluorescence microscope with a filter system for FITC and 400X magnification is required. The fluorescence pattern (form, density etc.) of the Negative and Positive Control is considered as reference pattern. Patterns of reactivity different than those seen in the controls must be considered non-specific, which means negative.

Positive test result at \geq 1:80 screening dilution: Leishmanie show a clear, yellow-green fluorescence on their membrane and flagella area. IgG titers of 1:80 or higher are considered to reflect infection at an

undetermined time. A further dilution of the positive samples is recommended to determine the highest dilution that is still positive (endpoint titer). An acute infection (2 to 4-fold titer increase) can only be determined by the titer determination of a coupled serum test (2 samples in an interval of 2-3 weeks). The interpretation of test results should always be based on anamnestic and especially clinical data and additional laboratory parameters.

Recommended cut-off 1:80:

Leishmanie (membrane, flagella) show a weak yellow-green fluorescence. The cut-off can vary depending on the region and the origin of the sample (dependent on the prevalence and on the state of endemicity). Therefore, it is recommended for every laboratory to determine an individual cut-off.

Negative test result at 1:80 screening dilution:

Leishmania do not show any yellow-green fluorescence, they are coloured greyish-red. Report as negative for Leishmania antibody.

Divergent fluorescent reaction:

Strong fluorescence of single (< 5 %) small, strongly globular Leishmania with shortened or missing flagellum among a majority of normally shaped non-fluorescent promastigotes must be considered as non-specific reaction and therefore as negative.

Possible cross reactions:


Due to immunological similarity to Trypanosoma cruzi, with sera from South and Central America, potentially false positive test results (antigen cross reactions) are possible.

QUALITY CONTROL

The negative control serum and the positive control serum should be assayed with each daily run. The negative control well is an example of a non-reactive serum, with either uniform red counterstain or slight, but uniform greenish staining. The fluorescence intensity of the positive control is an example of medium intensity fluorescence. If one of the controls does not react as specified, the assay run should be considered invalid, reagent components and procedural steps should be re-checked, and the assay repeated.

PRECAUTIONS

For veterinary use only. Follow the instructions carefully. Since no testing can assure the absence of infectious agents, these reagents, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly. Conjugate is photosensitive, for protection store in the dark and return to storage after use. Conjugate contains Evans blue dye, avoid ingestion and skin contact. Do not use after the expiration date.

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