Professional Veterinary Diagnostic Test
BioPronix product line
 For Veterinary use only
 For in vitro use only
 Manufactured by Agrolabo S.p.A.



INTERPRETATION OF THE RESULTS

The titre of IgG antibodies is determined by comparing the intensity of color between the Test line and Control line.

- If the T line (test line) has a lower intensity in comparison with C line (control line), the foal needs to receive an integration of immunoglobulins.
- If the T line (test line) has a less intensity in comparison with C line (control line), monitor and supplement the foal if it shows signs of disease or weakness.
- If the T line (test line) has the same intensity in comparison with C line (control line), the quantity of IgG is normal.
- If the T line (test line) has a higher intensity in comparison with C line (control line), the foal has received a very good transfer of IgG with maternal colostrum.

QUICK GUIDE

1. Transfer the filled loop by immersing the loop tip in the bottle with diluent
2. Mix sample and diluent by inverting the bottle (blue cap). With the pipette, dispense 1 drop of diluted whole blood, serum or plasma of foal in window marked S
3. Hold diluent bottle vertically and add 1 drop to the window marked S
4. Wait 5-7 minutes to read the results, not later than 10 minutes
5. Interpretation of the results

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
Rev. R1-27/04/22

TEST PRINCIPLE

FOAL IgG IC is an immunochromatographic test for the detection of immunoglobulins G (IgG) in foal whole blood, serum or plasma. The sample is absorbed into the cellulose pad. The IgG of foal move and bind with gold colloidal anti-IgG antibody present in the conjugate pad. The complex Ab-Ab will bind with the anti-IgG foal antibodies present on the membrane at level of line T. The result can be shown through the presence of a red line at the level of the T window (test line) and at the level of the C window (control line).

KIT COMPONENTS

- Disposable loop tips
- Disposable pipettes
- Vials for samples preparation
- Vial of diluent
- Foal IgG IC devices

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C). When properly stored, the product remains stable for its full shelf life. Keep the kit away from excessive heat/cold or direct sunlight.

PRECAUTIONS

IMPORTANT: open individual device packaging immediately prior to use. For veterinary use only. Do not use kit components after expiry date. The test device should be placed horizontally on a flat surface while testing. Use a new disposable pipette for each sample.

SAMPLES

Foal whole blood with anticoagulant, serum or plasma. All samples need to be brought to room temperature before use. A rapid identification of low IgG levels is very important for the early initiation of treatment of immunodeficient foals. IgG levels of foals should be measured when they are about 12-24 hours old.

TEST PROCEDURE

1. Remove the device from its individual wrapping and place it on a flat surface.
2. Insert the loop tip in the sample (whole blood, serum or plasma) and pick up the required volume of sample. Transfer the loop tip in the vial with diluent (blue cap). Mix the sample inverting the vial.
3. Holding the pipette vertically, dispense 1 drop of diluted sample into window S of the device.
4. Keeping the diluent bottle vertically, dispense 1 drop into window S of the device.
5. Read the results within 5-7 minutes and no later than 10 minutes.