**AGROLABO**  
Leishmania Ab ELISA 96 wells

# Leishmania Ab ELISA 96 wells

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# Leishmania Ab ELISA 96 wells

*Kit ELISA for the detection of specific antibodies to Leishmania infantum in canine serum or plasma samples*

## TEST PRINCIPLE

This test is performed as an indirect enzyme linked immunosorbent assay (indirect ELISA). The antigen of Leishmania is fixed on the wells of the plate. The samples to be analyzed are distributed in each well. If the sample contains specific antibodies anti-Leishmania, these antibodies will bind to the antigen adsorbed, forming the Ag-Ab complex. After washing to eliminate all non-fixed material, the presence of anti-Leishmania antibodies is detected by a specific anti-canine IgG monoclonal antibody, labelled with horseradish peroxidase (conjugate). All non-fixed material is removed by a second washing. A specific enzyme substrate is then added to the wells which causes a colour change reaction when the conjugate is present. The addition of stop solution no. 3 (fluoridric acid) allows a visual reading (blue colour), while the addition of stop solution no. 1 (sulphuric acid) allows a spectrophotometric reading (yellow colour).

## TEST COMPONENTS

- 1 x 96 wells plate (in strip of 8 wells);
- 1 dropper of positive control - ready for use;
- 1 dropper of negative control - ready for use;
- 1 dropper of conjugate - 100X;
- 1 bottle of washing solution - 10X;
- 1 bottle of diluent - ready for use;
- 1 dropper of substrate (TMB) - ready for use;
- 1 bottle containing stop solution no. 1 (sulphuric acid) - ready for use;
- 1 bottle containing stop solution no. 3 (fluoridric acid) - ready for use.

## MATERIALS REQUIRED BUT NOT SUPPLIED

Test tubes for sample dilutions; precision pipettes and tips; spectrophotometer at a wavelength of 450 nm.

## PRECAUTIONS

1. Read the instructions carefully.
2. Bring all reagents to room temperature (+20-25°C) prior to use.
3. Do not substitute or mix instructions or reagents from different kits or lots.
4. The substrate is very sensitive to the light and to contaminations. Do not use the tip directly in the vial. The stop solution is a diluted strong acid. Handle the two reagents with care.
5. Do not use components after expiry date.
6. Do not smoke, eat or drink where specimens or kit reagents are being handled.
7. In each utilization of the kit the positive and negative controls must be tested.
8. Use a new tip for each sample.
9. Observe each incubation time.
10. Prepare each sample according to instructions.

11. Avoid any contamination of the reagents.
12. All kit components must be stored at +2-8°C.

## WASHING STEPS

The washing steps may be carried out using an automated plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well. After the incubation, the washing steps should be carried out according to the following instructions:

- throw out the content of the plate by vigorously turning the plate over to avoid the possible mixture of the content from one well to another,
- dispense a volume of 300 µl of wash buffer into each well,
- shake the plate gently, avoiding contamination between wells,
- turn over the plate and tap the wells vigorously to remove the wash buffer,
- repeat the process as many times as is indicated by the kit's instructions,
- prior to emptying out the plate's contents during the final washing step, ensure that the next reagent to be added to the plate is ready to use,
- the wells should not be left dry for any longer than strictly necessary in order to dispense the reagent of the following step,
- after the last washing step, blot the empty wells face down on clean paper towels or absorbent paper to remove residual wash buffer.

## PREPARATION OF SAMPLES

Canine serum or plasma samples.

Prepare a **1:100** dilution of each sample to be tested using diluent.

## PREPARATION OF REAGENTS

### Wash buffer solution:

Dilute one part of the wash buffer supplied into 9 parts of distilled or deionized water (ie. 100 ml of concentrated wash buffer with 900 ml of water). Once diluted, the wash buffer solution may be stored at +2-8°C.

### Conjugate:

To be prepared immediately prior to use. Dilute the required amount of conjugate **1:100** with the diluent. Mix the solution prior to use. Prepare only the required amount and do not reuse surplus conjugate.

## TEST PROCEDURE

1. Bring all reagents to room temperature before use.
2. Add 100 µl of positive and negative controls in two different wells. Add 100 µl of diluted samples to the remainder wells of the plate. Delicately shake the plate.  
**Seal the plate and incubate at room temperature for 10 minutes.**
3. Wash cycle 4 times or more according to the procedure previously described.
4. Dispense 100 µl of previously prepared conjugate into each well. Delicately shake the plate.  
**Seal the plate and incubate at room temperature for 10 minutes.**
5. Repeat the wash cycle 4 times or more according to the procedure previously described.

6. Add 100 µl of substrate to each well. Delicately shake the plate.  
**Seal the plate and incubate at room temperature for 5 minutes.**
7. Add 100 µl of stop solution to each well.

**ATTENTION:** If you need a visual reading use stop solution no. 3 (fluoridric acid).  
If you need a spectrophotometric reading add the stop solution no. 1 (sulphuric acid) and read the absorbances of each well at 450 nm.

## READING AND INTERPRETATION OF THE RESULTS

The reading of the results differs depending on the stop solution used.

### VISUAL READING – Stop solution no. 3 (fluoridric acid)

**Positive samples:** a blue colour will develop. Colour intensity is proportional to the quantity of anti-Leishmania antibodies in the sample.

**Negative samples:** no colour will develop.

**Doubtful samples:** an intermediate colour light blue will indicate a low antibodies concentration; it is suggested to repeat the test after 3-4 weeks.

### SPECTROPHOTOMETRIC READING – Stop solution no. 1 (sulphuric acid)

The reading of the absorbance values (OD) must be carried out at a wavelength of 450 nm. If you run the samples or controls in duplicate wells, the OD value for the sample and controls will be the mean of both wells.

### Test Validation

The test is to be considered valid when:

OD of positive control is  $\geq 0.8$

OD of negative control < Negative Cut-off

### CUT-OFF value determination

Taking the positive value, will be determined the following cut-off:

Negative cut-off: OD positive control x 0.30

Positive cut-off: OD positive control x 0.35

### Interpretation of the results

**Positive Samples:** samples with an OD higher than the positive cut-off value.

**Negative Samples:** samples with an OD lower than the negative cut-off value.

**Doubtful Samples:** samples with an OD between both cut-offs. In these occasions, it is suggested to repeat the test after 3-4 weeks.

### TITRATION:

The antibody titre is represented by the highest dilution of the sample that has an absorbance value higher than the positive cut-off.

**AGROLABO**

**Manufactured by;**  
Agrolabo SpA

Head Office, Laboratories  
and Production Centre  
Diagnostic Division  
Via Masero 59  
10100 Scarmagno (TO)  
Italy  
Tel +39 0125 731111  
Fax +39 0125 731190  
E-mail [agrolabo@agrolabo.it](mailto:agrolabo@agrolabo.it)  
[www.agrolabo.it](http://www.agrolabo.it)  
[shop.agrolabo.it](http://shop.agrolabo.it)