Professional Veterinary Diagnostic Test
BioPronix product line
 For Veterinary use only
 For in vitro use only
 Manufactured by Agrolabo S.p.A.

INTERPRETATION OF THE RESULTS

Negative Result

The test is negative if a red line appears in window **C** (control line).

Positive Result

The test is positive if two red lines appear. One line band will appear in window **T** (test line) and one in window **C** (control line).

Invalid Result


The test is invalid if no coloured band appears in window **C** (control line), even if a coloured band appears in window **T** (test line).

QUICK GUIDE


1. Remove the device from its wrapping and place it on a flat surface. Use the pipette to add 1 drop of serum, plasma or whole blood with anticoagulant into window **S** of the device
2. Hold the diluent bottle in a vertical position and add 2 drops of diluent in window **S** of the device
3. Wait 10-15 minutes to read the results, not later than 15 minutes
4. Interpretation of the results

WARNING


Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

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
Leishmania IC - Rapid immunochromatographic test for the detection of antibodies to Leishmania

USER MANUAL

Rev. R4 -14/11/24

TEST PRINCIPLE

Leishmania IC is an immunochromatographic test for the detection of antibodies to *Leishmania infantum* in canine whole blood, serum or plasma.

If the sample contains specific antibodies against *Leishmania infantum* (positive sample), these antibodies will bind to the Leishmania antigen present on the membrane, forming the Ag-Ab complex at the level of the T window, developing a red line (test line).

The sample continues to migrate along the membrane up to window C, where a reagent able to bind the antibodies conjugated with colloidal gold is immobilized resulting in the formation of a red line (control line).

KIT COMPONENTS

- Disposable pipettes
- Vial with diluent
- LEISHMANIA IC devices

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C).

When properly stored, the product remains stable for its full shelf life.

Keep the kit away from excessive heat/cold or direct sunlight.

PRECAUTIONS

IMPORTANT: open individual device packaging immediately prior to use.

For veterinary use only.

Do not use kit components after expiry date.

The test device should be placed horizontally on a flat surface while testing.

Use a new disposable pipette for each sample.

SAMPLES

All samples need to be brought to room temperature before use. Previously frozen or stored samples must be centrifuged before use.

Use whole blood samples with anticoagulant (EDTA, heparin, etc), serum or plasma, fresh, or refrigerated (+2-8°C) for no more than 5 days. For longer storage, samples should be frozen.

Haemolyzed or lipemic samples may affect test results.

TEST PROCEDURE

1. Remove the device from its individual wrapping and place it on a flat surface.

2. Gently mix the sample by inverting. Holding the pipette vertically, add 1 drop of serum, plasma or whole blood with anticoagulant sample in the window **S** of the device. Wait for the sample to be absorbed.

3. Holding the diluent bottle vertically, dispense 2 drops of the diluent buffer in the window **S** of the device. **Note: In case of difficulty in sample migration, it is recommended to add 1 more drop of diluent in window S.**

4. Read the result after 10-15 minutes and no later than 15 minutes.